

FEATURED TECHNOLOGY

MANUALLY OPERATED DRIVER

SAM IO® is a manually operated intraosseous access system. Catheter placement is achieved by continuously actuating (repeatedly compressing) the driver's trigger assembly while gently guiding needle assembly into position. Repeated, full trigger actuation creates rotational spin of the needle assembly which, when combined with gentle downward pressure, results in controlled IO Needle placement.





ENGINEERED TO PRESERVE LIFE

SAM IO® NEEDLES

A stylet and catheter assembly secured by a press-fit connection. 15 mm, 25 mm, and 45 mm lengths are intended to address a broad range of patient demographics and multiple application sites. New visible stylet hub allows easy differentiation between stylet and catheter.

NEEDLEVISE BLOCK

Included in every Needle Set. Single-use sharps containment device for the stylet.

3 SAM IO® DRIVER

A light weight (2.5 oz) reusable manually operated intraosseous access driver made of a rugged Polycarbonate.

4 SAM IO® STABILIZER

Cylinder design intended to protect the catheter from dislodgement.

5 EXTENSION TUBING

Connects to the catheter for flexible infusion.

6 SAM IO®NEEDLE ADAPTOR

Attaches to the SAM IO® Driver to accommodate alternative IO needles.



PRODUCT SPECIFICATIONS

DRIVER SIZE & WEIGHT

Length: 5.5 in (14.0 cm) Width: 1 in (2.5 cm) Height: 4.25 in (10.8 cm) Weight: 2.5 oz (70.9 g)

* Size and weight vary by configuration and manufacturing process

SAM® is a registered trademark in the United States and other countries. $2025\,$ SAM Medical Products. All rights reserved.

MODELS AVAILABLE

Driver: 10700

NSN# 6515-01-694-6094

Stabilizer:

1-Pack IO740-1P NSN# 6515-01-694-6043 5-Pack IO740-5P

Adaptor:

10742

Needles:

15 mm 1-Pack IO705W-1P NSN# 6515-01-694-6069

15 mm 5-Pack IO705W-5P

NSN# 6515-01-694-6103 **25 mm** 1-Pack IO706W-1P

NSN# 6515-01-694-6098

25 mm 5-Pack IO706W-5P NSN# 6515-01-694-6052

45 mm 1-Pack IO707W-1P

NSN# 6515-01-694-6025 **45 mm** 5-Pack IO707W-5P

NSN# 6515-01-694-6112

GLOBAL SALES

+1.503.639.5474

sammedical.com

The CE marking on the SAM IO® Adaptor and SAM IO® Driver indicates self-declaration of conformity. No third-party certification or inspection has been conducted.







ZIEN Medical Technologies, Inc. 2490 South 300 West Salt Lake City, UT 84115 USA